NUK[®]

Thermometer 3 in 1

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Contents

1. Description	30
1.1 Intended use	30
1.2 Description of function	31
1.3 Scope of supply	31
1.4 Symbols on the device/packaging	31
2. Key to parts	33
3. Safety Instructions	33
3.1 Explanation of terms	33
3.2 General information	33
3.3 Special notes for this device	34
3.4 For the safety of your child	35
3.5 Handling batteries	35
3.6 Material damage	36
4. Starting the device	37
5. Operating the device	37
5.1 Basic information	37
5.2 Special things to consider when measuring a fever	38
5.3 Measurement in the ear	39
5.4 Measurement at the temple	40
5.5 Measuring surface temperatures	41
5.6 Measuring room temperature	43
5.7 Performing a follow-up measurement	43
5.8 Previous measurement	43
5.9 Memory (25 entries)	43
5.10 Switching between °C and °F	44
5.11 Silent mode	44
5.12 Switching off the device	44

6. Cleaning, maintaining and storing the device	45
6.1 Cleaning	45
6.2 Changing the battery	45
6.3 Storage	46
7. Troubleshooting	46
8. Technical data	48
8.1 Production date	50
9. Guidance and manufacturer's declaration	51
9.1 Electromagnetic emissions	51
9.2 Electromagnetic immunity	52
9.3 Recommended separation distances	55
10. Disposal	57
10.1 Equipment	57
10.2 Disposable/rechargeable batteries	57
10.3 Packaging	57
11. Warranty	58
11.1 Service Centres	59
11.2 Supplier	59
12. EU declaration of conformity	59

1. Description

1.1 Intended use

The infrared ear/temple thermometer NUK Thermometer 3 in 1 can measure the infrared heat generated by the temple or the

inner ear to reflect patient's body temperature accurately. The room temperature is shown on the display when it is idle.

Temperatures of surfaces can be simply and reliably measured.

This product should only be used in dry rooms. It must not be used out in the open or in humid rooms.

This product is intended for use in private households and is not suitable for commercial purposes.

1.2 Description of function

The thermometer uses an infrared sensor to measure the infrared radiation at the eardrum membrane or at the temple.

This radiation is converted into body temperature as a function of measurement either in the ear or at the temple, and is shown in the display.

The practical fever indication indicates an increased temperature using a red LED signal.

The memory makes it possible to save measurements and observe temperatures over time.

With its integrated battery, the thermometer is immediately ready for use.

1.3 Scope of supply

- 1 thermometer, consisting of:
 - 1 thermometer
 - 1 cover
 - 1 battery CR2032 (inserted into device)
 - 1 stand
- 1 set of operating instructions

1.4 Symbols on the device/ packaging

Symbol	Meaning
	Please read the instructions for use in full before using the device for the first time.
	Caution
IP22	Protection level for water and dust

Symbol	Meaning
	Symbol for permitted temperature range when storing the device
	Paper Recycling
	Authorized representative in the European community
	Manufacturer
	Medical device of type BF
	The CE mark and Notified Body Registration Numbers, the requirement of Annex II from Medical Device Directive 93/42/EEC are met.
	Battery Recycling

Symbol	Meaning
	Indicates this device is subject to the Waste Electrical and Electronic Equipment Directive in the European Union. To protect the environment, dispose of useless device at appropriate collection sites according to national or local regulations.

2. Key to parts

- | | | | |
|---|-------------|----|---------------------------------------|
| 1 | Cover | 6 | Unlocking opening |
| 2 | SCAN button | 7 | Stand-by (ON/MEM) button |
| 3 | LCD display | 8 | Battery CR2032 (inserted into device) |
| 4 | Sensor lens | 9 | Battery compartment lid |
| 5 | Probe head | 10 | Stand |

3. Safety Instructions

3.1 Explanation of terms

The following keywords are used in these user instructions:

Warning!

High risk: ignoring this warning could result in a risk to life and limb.

Caution!

Medium risk: ignoring this warning may result in damage.

Note:

Low risk: points to be noted when handling the device.

3.2 General information

- Please read these instructions for use carefully before using the device. These instructions form part of the device and must be kept to hand at all times.
- Use the device and its accessories only for the proper use described (see “1.1 Intended use”).
- Do not use any accessory parts that are not expressly approved for the device.
- Always check for any visible damage before using the device. If you discover any damage, do not use the device.

- Changes may only be made to the device with the consent of the manufacturer.
- Always seek the advice of a doctor in the event of high fever or ambiguous measurements.

3.3 Special notes for this device

- This device is not designed to be used by people (including children) with physical, sensory or mental disabilities or with insufficient experience and/or knowledge, unless a person responsible for their safety is supervising them or instructing them on how to use the device.
- Do not use any additional probe covers on the device.
- If the probe head on the device glows red when the temperature is measured, continue to monitor the body temperature constantly. If in doubt about the person's condition, seek medical assistance.
- Clean the device after every use in order to prevent any reinfection or infection of other persons.
- The "Surface temperature" setting cannot be used to measure body temperatures/fevers.
- This device can be used by children aged 11 and up if they have at least 5 years' intensive experience of reading and if they are under supervision or have been briefed on the use of the device.

- Check the function and accuracy of the device after any heavy impact (e.g. falling from a height of more than 1 metre) by taking several measurements. If you have any doubts concerning the measurement results, please contact your retailer.
- The use of the device doesn't replace the investigations of the doctor.
- In order to avoid inaccurate measuring, the probe head **5** and the sensor lens **4** needs to be cleaned **after every use** (see "6.1 Cleaning" on page 45).

3.4 For the safety of your child

Warning! Children are often unaware of danger and are injured as a result. Please note therefore:

- Only allow children to use this device under the supervision of an adult.
- Always ensure that the device is kept out of reach of children when not in use.
- Ensure that the packing foil does not become a deadly trap for a child. Packing foils are not toys.

3.5 Handling batteries

Warning! To prevent damage to health:

- To prevent a risk of explosion, normal batteries must not be charged, heated or disposed of by burning.
- Never try to open batteries.
- To prevent injury, if battery acid has leaked avoid contact with skin, eyes and

mucous membranes. In the event of contact with the acid, rinse the affected areas immediately with plenty of clean water and contact a doctor at once.

- Keep the battery out of reach of children. There is risk of ingestion. If there is a suspicion of swallowing the battery, immediately see a doctor.
- Batteries may contain toxic substances, so follow the instructions for disposal (see “10.2 Disposable/rechargeable batteries“ on page 57).

Caution! To prevent material damage:

- The batteries must be removed when flat or if the device is not going to be used for an extended period, otherwise leaking batteries may damage the device.
- Be sure only to use type CR2032 batteries.

3.6 Material damage

Caution! To prevent material damage, follow the instructions below.

- Do not expose the device to direct sunlight or extreme heat, as UV rays and overheating may make the plastics brittle and the electronics might be damaged.
- Never use the device in very damp or wet surroundings.
- Never use caustic or abrasive cleaning agents, as this could damage your device.
- Do not submerge the device in water or expose it to direct moisture.

4. Starting the device

1. Remove all packaging material from the device.
2. Pull out the insulating strip of the battery.
The thermometer switches on with a beep.

Note: If the insulating strip tears off when you pull it out, take out the battery and remove the rest of the insulating strip (see “6.2 Changing the battery” on page 45).

3. Clean the device (see “6.1 Cleaning” on page 45).

5. Operating the device

5.1 Basic information

Fever is a state of increased body temperature.

Body temperature:

A “normal” temperature is around 36.5 to 37.0 °C (97.7 – 98.6 °F). The highest temperature will be read in the early evening. The values differ according to where the temperature was taken (mouth, armpit, anus, ear, temple).

Causes:

Fever may result from colds, other illnesses

such as respiratory infections, some intestinal infections, hyperthyroidism or the presence of tumours.

Symptoms and diagnoses:

Sweating, shivering or chills are indications of an increased body temperature.

You should contact a doctor:

When fever starts, it may have a negative effect on the body, depending on age, the constitution and the current well-being of the person. The use of the device doesn't replace the investigations of the doctor.

Please inform about the critical values and consult a doctor if these are exceeded.

5.2 Special things to consider when measuring a fever

Please note that when taking the temperature several times in succession from the same ear, the results may fluctuate by up to half a degree. These results may fluctuate to an even greater extent should you change to the other ear.

In the event of an ear infection, damage to the ear drum or too much earwax, take the temperature at the temple instead.

If the person has been lying on their ear for a longer period of time, the other ear should be used.

The ear canal of infants less than six months old is often very narrow, which may lead to values that are too low. It is therefore recommended that you take the temperature of infants who are less than six months old at the temple.

Factors which may affect the measured value when **measuring in the ear**:

- The bend of the ear canal
- Ear wax
- Incorrect positioning of the thermometer at the ear drum
- Moisture in the ear canal

Factors which may affect the measured value when **measuring at the temple**:

- Thickness of the skin at the temple
- Sweat on the temple
- Vasoconstrictive medication
- Skin complaints
- Moisturising cream on the skin

The device must stay in stable ambient (room) temperature for 15 minutes before operating.

Do not hold the thermometer in your hands for longer than necessary. Otherwise, the device


will warm up and may read a temperature that is too low when checking for fever.

The device is intended for indoor use, where possible, try to avoid using the thermometer outdoors. External influences, such as the wind and outdoor temperature, may have an impact on the measurement.

If the temperature taken by the thermometer exceeds or equals 37.5 °C (99.5 °F), a red light will appear in the probe head to warn

the user of a possible fever. If the temperature taken is below 37.5 °C (99.5 °F), a green light will appear.

The person who is having their temperature taken should avoid any fluctuating ambient temperatures for at least 5 minutes before the measurement. Also avoid any strenuous physical activity or taking a bath for at least 30 minutes before the measurement.

The temperature of a healthy person may differ by 0.2 - 1.0 °C / 0.4 - 1.8 °F between the temple and the ear

5.3 Measurement in the ear

Note: The mode "Taking the temperature at the ear" is the standard mode after switching on the device.

Take the temperature a total of three times from the same ear. If the results differ, take for safety the highest temperature value.

1. Press the Standby (ON/MEM) button **7** to switch on the device.

The adjacent symbol will appear in the LCD display **3** and you will hear two beeps.

2. Remove the cover **1** from the probe head **5** (see figure **C**).
3. Gently pull back the ear to lessen the bend of the ear canal (see figure **D**).

4. Carefully guide the probe head **5** into the ear canal until it rests on the ear canal (see figure **E** and **F**).
5. Press and hold the SCAN button **2** until you hear a beep.
6. Remove the probe head **5** from the ear canal.

The egg-timer symbol is displayed and stops flashing after about 4 seconds. Two short beeps indicate that the result has been determined.

7. Read off the measured result in the LCD display **3**.

The device automatically switches to room temperature mode after approx. one minute.

5.4 Measurement at the temple

Make sure the temples are dry, clean and free from make-up when taking the

temperature. Avoid placing the device on any scarred areas.

1. Press the Standby (ON/MEM) button **7** to switch on the device.

The adjacent symbol will appear in the LCD display **3** and you will hear two beeps.

2. Remove the cover **1** from the probe head **5** (see figure **C**).
3. Press and hold the Stand-by (ON/MEM) button **7** and, at the same time, press the SCAN button **2** once.

The temples symbol will appear on the LCD display **3**. The "temple measurement mode" is thus activated, used for reading the temperature at the temples.

4. Touch the skin on the right or left temple with the sensor lens **4** of the thermometer (see figure **G**).

5. Briefly press the SCAN button **2** and gently move the probe head **5** in a circular motion over the skin on the temple (see figure **H**).

As you are doing this, you will hear a sound that indicates the temperature is being read.

Note: The time required to take the temperature ranges from 5 to a maximum of 8 seconds. This depends on how long the device needs to determine the correct temperature of the temples.


As soon as you hear a second sound, the measuring process is complete. The egg-timer symbol  is displayed and stops flashing after about 4 seconds. Two short beeps indicate that the result has been determined.

6. The temperature taken is shown on the LCD display **3**.

The device automatically switches to room temperature mode after approx. one minute.

5.5 Measuring surface temperatures

The thermometer has a measuring range of -22.0 to 80.0 °C (-7.6 to 176.0 °F) for measuring the temperature of

surfaces. This enables you to measure the surface temperature of radiators, cooker hotplates, ovens or similar surfaces, as well as liquids.

Note: Please note that steaming surfaces, such as water, food or milk, for example, may cause the sensor lens **4** to mist up, resulting in a measurement which may be inaccurate.

Warning! The “Surface temperature” setting cannot be used to measure body temperatures/fevers.

Note: Please note when measuring surfaces that it really is only the temperature of the surface which is determined, not the temperature of what is inside. There may be differences in temperature between the surface and the core, especially when measuring liquids. It may be much higher. Incorrect conclusions may lead to burns or scalds. Before measuring a liquid, stir it to even out these differences as far as possible and then take the measurement immediately.

1. Press the Standby (ON/MEM) button **7** to switch on the device.

2. Hold down the Standby (ON/MEM) button **7** and press the SCAN button **2** at the same time until the adjacent symbol  appears in the LCD display **3**.

This activates “Surface measurement” mode which enables you to determine the temperature of surfaces.

3. Remove the cover **1** from the probe head **5** (see Figure **C**).

4. Hold the sensor lens **4** of the thermometer as close as possible to the surface to be measured without contacting it directly.

Note: Under no circumstances may the thermometer be immersed in liquids or contact hot cooker hotplates.

5. Briefly press the SCAN button **2** to perform the measurement.

6. Read off the temperature measured in the LCD display **3**.

7. Press and hold the SCAN button **2** to update the measurement constantly.

The device automatically switches to “Room temperature” mode after approx. one minute.

5.6 Measuring room temperature

If the thermometer is not operated for around 1 minute, the device automatically switches to “Room temperature” mode.

The symbol (🔧) and the room temperature measured appear on the LCD display **3**.

5.7 Performing a follow-up measurement

After you have performed one measurement, you can perform a follow-up measurement by keeping the SCAN button **2** depressed again.

5.8 Previous measurement

For every new temperature measurement, the result from the previous measurement will be shown for comparison at the top of the LCD display **3**.

The temperature measured previously is displayed regardless of whether it was measured in the ear or at the temple.

5.9 Memory (25 entries)

1. With the device switched on, press the Standby (ON/MEM) button **7** simultaneously.

Memory 1 appears and displays the last temperature measurement when the button is released.

2. Pressing the button again in this way displays the previous 25 temperature measurements performed.

The memory is represented in the LCD display **3** by the adjacent symbol. 

The appropriate symbols also indicate whether temperature was measured in the ear or at the temple. 

5.10 Switching between °C and °F

1. If the device is in "Room temperature" mode, press and hold the SCAN button **2** whilst simultaneously pressing the ON/MEM button **7** for 3 seconds.

The character in the LCD display **3** changes: "°C" to "°F".

2. The same procedure applies when switching from °F to °C.

5.11 Silent mode

You can switch off the sounds using silent mode.

1. With the device switched on, press Standby (ON/MEM) button **7**.

The silent mode symbol starts flashing in the LCD display **3**. 

Release Standby (ON/MEM) button **7** again to activate silent mode.

The same procedure applies when deactivating silent mode.

5.12 Switching off the device

To switch off the thermometer, press and hold the Standby (ON/MEM) button **7** until "OFF" appears in the LCD display **3**.

6. Cleaning, maintaining and storing the device

6.1 Cleaning

Caution!

- In order to avoid inaccurate measuring, the probe head and the sensor lens needs to be cleaned **after every use**.
 - Do not submerge the device in water or expose it to direct moisture.
 - The probe head **5** is the most sensitive part of the thermometer. Treat this part of the device with particular care when cleaning the sensor lens **4**.
1. For cleaning **only use** a cotton bud dipped in a solution of app. 70% alcohol (Ethanol) (see figure I).


Caution!

- Do not use other chemical agents or disinfectants to clean as this may damage the device.

- Make sure here that the tip of the device is pointing downwards so that none of the cleaning fluid enters the housing.
2. Then leave the probe head to dry completely for at least one minute.

6.2 Changing the battery

The battery must be replaced as soon as the following symbol appears in the LCD display **3**.

1. Take a small, thin object (e.g. a paper clip) and press it into the unlocking opening **6** of the battery compartment lid **9** while, at the same time, pushing this opening outwards with both thumbs (see figure **K**).
2. Firmly holding the device, carefully lift out the battery **8** using a small, thin object (e.g. a paper clip) (see figure **L**).

3. Insert the new CR2032 battery with the positive terminal facing upwards and the minus terminal facing downwards by sliding the battery underneath the metal hook and pressing down until it clicks into place (see figure **M**).

Caution! If incorrectly fitted, the battery compartment and even the device may become damaged.

4. Push the battery compartment lid **9** back onto the device until it audibly clicks into place.

6.3 Storage


We recommend that you retain the packaging of the thermometer so that the device can be packaged securely for transport.

If the device is not to be used for an extended period, we recommend switching off the device (see “5.12 Switching off the device“ on page 44).

Only operate and store the device in accordance with the conditions specified. (see “8. Technical data“ on page 48).

7. Troubleshooting

Error message	Problem	Possible cause	Solution
Er	The device is not working properly.	Software error	Remove the battery, wait one minute and then replace the battery (see sect. 6.2). If the error message appears again, contact your retailer.

	While the symbol  is flashing, a new measurement cannot be started in the ear/temple mode.	The device is not yet ready to take a new measurement	Wait until the symbol has stopped flashing. When it has stopped flashing, then start the new measurement.
AbH	Room temperature greater than 40.0 °C (104.0 °F).	The device was not stored as instructed	Place the thermometer in a room with a temperature of between 10.0 °C (50.0 °F) and 40.0 °C (104.0 °F) for at least 30 minutes.
AbL	Room temperature less than 10.0 °C (50.0 °F).		
	The device will not switch on and display anything.	Insulation strips have not been removed from the battery	Remove the insulation strips from the battery (see sect. 4.).
		Battery empty	Replace the battery (see sect. 6.2).

Hi	<p>(1) <i>In ear/temple mode:</i> Measurement temperature > +42.2 °C (108.0 °F)</p> <p>(2) <i>Surface temperature:</i> Measurement temperature > +80.0 °C (176.0 °F)</p>	Object to be tested outside the specified measurement range (too hot)	Take a measurement within the specified temperature range. In the event of a malfunction: Contact your retailer.
Lo	<p>(1) <i>In ear/temple mode:</i> Measurement temperature < +34.0 °C (93.2 °F)</p> <p>(2) <i>Surface temperature:</i> Measurement temperature < -22.0 °C (-7.6 °F)</p>	Object to be tested outside the specified measurement range (too cold)	

8. Technical data

Temperature range:	
Ear and temple mode:	+34,0~42,2 °C (93,2~108,0 °F)
Surface temperature mode:	-22,0~80,0 °C (-7,6~176,0 °F)
Operating temperature:	+10,0~40,0 °C (50,0~104,0 °F), 15%~85% RH (non-condensing)
Storage temperature:	-20,0~+50,0 °C (-4,0~122,0 °F), relative humidity ≤85%
Transport temperature:	<70,0 °C (158,0 °F), relative humidity ≤95%

Accuracy:	
Ear mode:	$\pm 0.2\text{ }^{\circ}\text{C} / 0.4\text{ }^{\circ}\text{F}$ in the measuring range $35.5\sim 42.0\text{ }^{\circ}\text{C}$ $(95.9\sim 107.6\text{ }^{\circ}\text{F})$ $\pm 0.3\text{ }^{\circ}\text{C} / 0.5\text{ }^{\circ}\text{F}$ outside this range Clinical Bias: $-0.5 \sim -0.7\text{ }^{\circ}\text{C}$ Limits of Agreement: 1.02 Repeatability: $0.19\text{ }^{\circ}\text{C}$
Temple mode:	$\pm 0.2\text{ }^{\circ}\text{C} / 0.4\text{ }^{\circ}\text{F}$ in the measuring range $35.5\sim 42.0\text{ }^{\circ}\text{C}$ $(95.9\sim 107.6\text{ }^{\circ}\text{F})$ $\pm 0.3\text{ }^{\circ}\text{C} / 0.5\text{ }^{\circ}\text{F}$ outside this range Clinical Bias: $-1.3 \sim -1.5\text{ }^{\circ}\text{C}$ Limits of Agreement: 1.47 Repeatability: $0.22\text{ }^{\circ}\text{C}$
Surface temperature mode:	$\pm 0.3\text{ }^{\circ}\text{C} / 0.5\text{ }^{\circ}\text{F}$ in range $22.0\sim 42.2\text{ }^{\circ}\text{C}$ ($71.6\sim 108.0\text{ }^{\circ}\text{F}$), otherwise $\pm 4\%$ (however at least $\pm 2.0\text{ }^{\circ}\text{C} / 4.0\text{ }^{\circ}\text{F}$)
Atmospheric pressure:	800~1013 hPa
Protection class:	BF
Degree of protection:	IP22

Battery:	1 CR2032 lithium battery
Battery life:	approx. 3,000 measurements (~1 year)
Comply with:	ASTM E1965-98, EN ISO 80601-2-56, IEC/EN60601-1-2(EMC), IEC/EN60601-1(Safety standards), ISO10993, RoHS.

Note: The thermometer converts the ear/temple temperature into its “oral equivalent” (referring to the results of the clinical assessment).

8.1 Production date

The production date is contained in the serial number:

Example: Ex. SN: E512A000001

Ex. SN: “Serial number”

E: code for “External”

5: last digit of production year, in this case, 2015

12: production month, in this case, December

A000001: consecutive serial number

9. Guidance and manufacturer's declaration

9.1 Electromagnetic emissions

The THixyz series is intended for use in the electromagnetic environment specified below. The customer or the user of the THixyz series should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The THixyz series uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The THixyz series is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

9.2 Electromagnetic immunity

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of the THixyz series, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	<p>Recommended separation distance</p> $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2,3 \sqrt{P} \text{ 800 MHz to 2,5 GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the THxyz series is used exceeds the applicable RF compliance level above, the THxyz series should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the THxyz series.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	6 kV contact 8 kV air	6 kV contact 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Surge IEC 61000-4-5	1 kV line(s) to line(s) 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
interruptions and voltage variations on power suppl	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the THixyz series requires continued operation during power mains interruptions, it is recommended that the THixyz series be powered from an uninterruptible power supply or a battery.
input lines IEC 61000-4-11	70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec		

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: UT is the a.c. mains voltage prior to application of the test level.

9.3 Recommended separation distances

between portable and mobile RF communications equipment and the ME equipment or ME system.

The THixyz series is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the THixyz series can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the THixyz series as

recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter in W	Separation distance according to frequency of transmitter in m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	N/A	0,12	0,23
0,1	N/A	0,38	0,73
1	N/A	1,2	2,3
10	N/A	3,8	7,3
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

10. Disposal

10.1 Equipment

The symbol of the crossed-through wheellie bin means, within the European Union, that the product is subject to a separate waste collection. This applies to the product and all accessories marked with this symbol. Products labelled as such must not be disposed of with the normal household waste collections; they must be taken to a recycling bank for electrical and electronic equipment. Recycling helps to reduce the use of raw materials and is eco-friendly.

10.2 Disposable/rechargeable batteries

Batteries must not be disposed of with domestic waste. Flat batteries must be disposed of properly. Special containers for battery disposal are provided for this purpose at shops selling batteries or at municipal collection points.

Batteries and rechargeable batteries on which the following letters are printed contain, among other things, the following harmful substances: Cd (cadmium), Hg (mercury), Pb (lead).

10.3 Packaging

When you wish to dispose of the packaging, please observe the relevant environmental regulations for your country.

11. Warranty

In addition to the seller's warranty obligations arising from the contract to purchase, we as the manufacturers provide a warranty of 24 months from purchase of the device as long as the device is used properly and the operating instructions are followed. Evidence of the date of purchase and the type of device should be provided by means of a receipt.

The purchaser's statutory rights are unaffected by this manufacturer warranty.

During the warranty period, we undertake to remedy any defects attributable to poor materials or workmanship. Worn parts and batteries are excluded from the warranty.

Minor deviations from specified quality which do not affect the value or usability of the device are not covered by the warranty. We are likewise unable to accept the warranty if the defects on the device are due to transport damage for which we are

not responsible, if they are attributable to misuse or lack of care or if the device has been interfered with by persons not authorized by us to do so.

The warranty service will involve repair, replacement of parts or replacement of the device at our discretion. The warranty period will be neither extended nor restarted by the performance of warranty services. The warranty for any spare parts fitted will expire with the warranty for the device as a whole.

Extended or different claims, especially those relating to the making good of damage sustained outside the device, are excluded unless liability is compulsory in law.

We will not accept the costs or risks of transport.

An device sent in without evidence of the purchase date will be treated as a repair.

The device will only be repaired following discussion with the customer.

Retain address, article number and LOT No. for possible future reference. The lot and serial number of the device can be found on the inside of the battery cover.

11.1 Service Centres

The addresses for service can be found under “NUK Contact addresses” on page 114.

Article no.: 10.256.345

11.2 Supplier

Please note that the address below **is not a service address**. In the event of problems or queries relating to your product, contact one of the Service Centres (see “NUK Contact addresses” on page 114).

MAPA GmbH
Industriestraße 21-25
27404 Zeven
Germany
www.nuk.com

12. EU declaration of conformity

The device meets all applicable European directives, as well as their respective essential requirements. These are detailed in the EU declaration of conformity, which can be requested from the manufacturer. The declaration of conformity can be found at www.nuk.de.

The thermometer meets the requirements of the standard DIN EN 12470-5 “Clinical thermometers, Part 5: Performance of infra-red ear thermometers (with maximum device).”

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NUK is a registered trademark of
NUK is a registered trademark of MAPA GmbH, Germany
MAPA GmbH, Germany