NUK®

Beckenbodentrainer Pelvic Floor Trainer

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Operating Instructions

D **GB** **F** **I** **E** **P** **PL** **RUS** **GR** **TR**

Rev. 400
Art.-Nr. 10.256.366

MAPA GmbH
Industriest. 21-25
27404 Zeven
Germany

www.nuk.com

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MAPA GmbH, Industriestraße 21-25, 27404 Zeven, Germany, www.nuk.com

(GB) Instructions for use

By buying the NUK pelvic floor trainer, you have acquired a product which will support postnatal recovery and contribute to strengthening the muscles of your pelvic floor. Please read these instructions through carefully before using for the first time.

Function

The pelvic floor is put under considerable strain during pregnancy and birth, and this can result in it no longer performing its function in full. This pleasantly small and silent training aid stimulates, strengthens and trains the muscles of the pelvic floor. It contains a lightweight ball which vibrates along with every movement you make; it is this stimulation which supports active muscle build-up. No more than a few minutes' training twice a day are required.

Cleaning

Before using for the first time, and after every use, the NUK pelvic floor trainer should be washed using only warm water and a mild soap. If you would like to disinfect the trainer, we recommend NUK Hygiene Spray.

Use

The NUK pelvic floor trainer should only be used after the birth once any birth-related injuries have healed, the lochia has stopped and postnatal exercising has begun.

Please check the pelvic floor trainer for visible signs of damage before using it each time. The NUK pelvic floor trainer is inserted lengthwise into the vagina, held by the end with the retrieval string, in a similar way to a tampon. For easier insertion and removal, please use a water-based lubricant - like the one from Billy Boy, for example. Do not use greasy or oily products. The elastic retrieval string remains visible and within reach outside the body. Depending on how you feel, the trainer can remain in your body for just a few minutes or for several hours. For best results, it should be worn during everyday physical activities. To facilitate removal, relax your pelvic floor and pull the NUK pelvic floor trainer out carefully with the aid of the retrieval string.

Keep the NUK pelvic floor trainer in a clean place. Do not expose it to direct sunlight.

Please note that the NUK pelvic floor trainer should not be used during your period, as the tissue is highly sensitive during this time. Women who would like to exercise their pelvic floor muscles during pregnancy should only use the NUK pelvic floor trainer after consulting their midwife or obstetrician.

Warning

The NUK pelvic floor trainer is not intended for anal use, as it may be retained by the anus. Please remove before vaginal intercourse.

The NUK pelvic floor trainer is not a toy and should be kept out of the reach of children.

Avoid damaging the trainer with fingernails, jewellery or other sharp objects.

In the event of vaginal infections or if the product is damaged, please dispose of the NUK pelvic floor trainer for safety reasons.

Retain the address and article no. in case of a complaint.

Distributed in the UK by: MAPA Spontex UK Ltd., Berkeley Business Park, Wainwright Road, Worcester WR4 9ZS, Great Britain, www.nukbaby.co.uk
Distributed in Ireland by: Intrapharma Consumer Health, Magna Business Park, Citywest Road, Dublin 24, Ireland

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NUK is a registered trademark of